

Frequency Synthesizer F-SCAN² SATELLITE

INSTRUCTION MANUAL



WARNING: User's with PACE MAKERS or PACE MAKER ELECTRODES should consult a cardiologist prior to use of this device. The device could perhaps interfere with – or even damage – the PACE MAKER.

WARNING: The device generates frequencies. It must be used as directed by the therapist. The use of cables, adapters or accessories other than the ones supplied damages the device and may cause burns or injuries.

WARNING: Never apply electrodes on open wounds.

Use of the SATELLITE

- Connect the application cable with its plugs marked red and blue to the electrodes or other application parts assigned by the therapist.
- Connect the application cable with its plug marked black to the socket OUTPUT of the SATELLITE.
- Set the program-switch (SELECT) of the SATELLITE to the position (between 1 and 10) of the application to run.
- Switch the SATELLITE ON (Press the button POWER and hold it down until a clicking sound is heard after which the button must be released).
- The device issues 20 clicking sounds (takes about 15 seconds) and starts the application without a pause. The electrodes or other application parts must be positioned during that time. The green LED above the button POWER is lit permanently while the clicking sounds are issued and flashes when a frequency is active.

The device signals acoustically the changeover from one frequency to the next and switches OFF at the end of the application. The green LED above the OUTPUT is lit when a frequency is active.